

D0A-144		Multi drugtest 4	
Abbreviation	Cut-Off mg/mL	Substances	Approximate detection times
MOP	300	pholcodine, heroin, codeine derivatives, levofanol, meperidine, morphine derivatives, nalofin	1-5 days
BZO	300	antham, acepam, acepates, azolam, diazepoxides	2-14 days
AMP	1000	amphetamine, phenylethylamine, phentermine, MDA, tryptamine, tyramine	1-6 days
THC	25	Cannabis (hashish, marijuana)	5-50 days after their use depending on the amount and regularity of use